



Ultraview SL3800

Spacelabs Healthcare - Connecting Innovation with Care

Spacelabs Healthcare is fully committed to making strong connections: between clinicians and their patients; between technological innovations and real-life patient care; between clinical needs and cost-efficient solutions. We are a dynamic company with decades of experience providing innovative products and services to healthcare professionals.

Our goals are simple:

Keep it easy to use, focus on system flexibility and scalability, and provide access to information anytime, anywhere. The result will continue to fulfill your greatest need: more time to care for your patients.



Scalable Central Surveillance

The Ultraview SL™3800 combines central monitoring for hardwired and wireless patient monitors, as well as patient-worn telemetry systems. Interactive network controls allow you to view, control, review and record information for any monitored patient, while providing alarm notification for all patients automatically. Scalable solutions like WinDNA® and the Intesys® Clinical Suite (ICS) save you valuable time by centralizing access to and management of patient information.

Standard features include: Graphic and Tabular Trends and interactive networking. Basic configuration provides eight waveforms; twelve and sixteen waveforms are optional. Multiple flat panel display sizes (15" or 19") are available.



Advantages

Consistency and flexibility for superior ease of use and optimal performance in any care environment.

- Default alarm limits and customizable parameter settings
- Modular, scalable architecture
- Compatible with Ultraview®, Ultraview SL and PCMS™ monitors and modules
- Capability for configuring as a bedside monitor
- Easily upgraded and/or maintained in the field

Integrated alarm management solutions to streamline event notifications and documentation.

- Alarm Limit Review provides a printable summary of alarm limits for all active parameters
- Flexport® system interfaces link alarms from standalone devices at the bedside
- Tone Access controls
- Alarm Watch delivers bed-to-central alarm notification




Superior network solutions deliver instant access to patient information.

- Centrally view, control, review and record patient data for all monitored patients
- WinDNA for delivering charting, lab, Internet, intranet and HIS applications
- Intesys Clinical Suite (ICS) provides departmental or enterprise solutions for Full Disclosure, Electronic Flowsheet and network laser printing
- HL7 interfaces to HIS applications

Ultraview SL3800



Collecting and centralizing information provides monitoring flexibility for every department.



Connecting Innovation with Care

Network

The Spacelabs monitoring network supports seamless data acquisition and exchange across the enterprise, addressing the clinician's need for continuous information management.

Supplies and Accessories

High-quality supplies and accessories are important for monitoring system performance and for economical and efficient patient care. All Spacelabs supplies and accessories are carefully designed and tested to guarantee optimal system performance.

Features

- Supports up to 16 waveforms and 26 parameters
- Flexible display options
- Compatible with all Ultraview and PCMS/PCIS monitors and telemetry systems
- Field upgradeable/repairable
- Reconfigurable as a bedside monitor

Options

- Drug Dose Calculations
- Vital Sign Calculations
- Full Bed Review (FBR)
- WinDNA

Biomed Spare Options

(for use in bedside mode)

- Data Shuttle®
- Perioperative Software
- Full View ECG Display
- Patient Data Logger™ (PDL)

Available accessories include:

- Hardwired or wireless keyboard and mouse
- USB adapter cable for PS/2 accessories
- GCX mounting and mobility solutions
- Intesys Clinical Suite (ICS)
- Flat panel displays: 15", 19"
- 2/4 channel recorder