

Flexports 90XXX and 91XXX

05-January-2018

077-0384-00 Rev A

CUSTOMER SERVICE NOTES

Product Covered: Flexports, 90XXX and 91XXX

Purpose: Third Party Software Versions Validated with Spacelabs Flexport Products

Description:

The Spacelabs Flexport products allow presentation of non-Spacelabs medical device data as an additional parameter on your Spacelabs monitoring system. Each non-Spacelabs medical device software version uses a unique communication protocol to interact with your Spacelabs medical device. The purpose of the protocol is to present data in a specific predictable format.

The list below identifies third party software versions validated for use with the Spacelabs Flexport products. Any other third party software version may not present data appropriately.

Manufacturer	Model	Firmware Version	Flexport Model	Spacelabs Firmware Version
BD/CareFusion/Viasys	Vela	2.02.08	91436C	1.00.02-002
BD/CareFusion/Viasys	Avea	3.9B	91436C	1.00.02-002
Drager	Apollo/Primus/Primus IE	V.4.50 (or4.51) Medibus 6.00 Medibus.X Enable	91443C	1.00.01
Drager	Perseus	V.1.12 Operating System, Medibus.X Protocol	91443C	1.00.01
Edwards	Vigilance II	V.0056	91439B	1.00.01-002
Edwards	Vigileo	V.01.10	91439B	1.00.01-002
Edwards	EV1000	V.01.10	91439B	1.00.01-002
GE Healthcare	Engstrom Carestation	5.05	91436E	1.00.01-004
GE Healthcare	Avance 5.0	5.0	91436D	1.00.02-002
Hamilton Medical	Galileo	3.2	90436A-07	1.00.02
Hamilton Medical	Raphael	3.2	90436A-07	1.00.02
Maquet	Servo-I	1,2,3,4	91436B	1.00.05-001
Masimo	Radical	Masimo Set-V4.1.0.1, Handheld-R4.2.0.1	90442A-02	1.00.52
		D-station (Radnet) R5.1.3.2	90436A-15	2.00.02EN

Manufacturer	Model	Firmware Version	Flexport Model	Spacelabs Firmware Version
Medtronic/Covidien/Puritan Bennett	PB840	4-070212 85-N	90436A-15	1.00.06EN
Oridion	Vital Cap Microstream	2.75	90442A-02	1.00.00EN
Pulsion	<i>Picco I/Picco Plus</i>	D7.0.7 (Demo Mode)	90439	1.00.01
Pulsion	Picco	3.1.0.8A	91439A	1.00.00-013
SenTec	Sentec	7.02	90442A-03	2.00.02EN
Spacelabs	Sirius	700900 9.62 Control Board 70090010.06 FrontPanel	90519B	1.00.05
Spacelabs	Blease	700900 9.62 Control Board 700900 9.02 FrontPanel	90519B	1.00.05
Spacelabs	Arkon	Contains a protocol Check	91443A	1.00.00-016

Contact your local Sales Representative for any further questions.