

# SUPPLIES & ACCESSORIES

## CARDIAC OUTPUT










# CARDIAC OUTPUT

## Cardiac Output Supplies

A variety of supplies are available through Spacelabs for 90496/91496 modules configured for cardiac output. These supplies facilitate the thermodilution method. Spacelabs offers adapters, cables, reference and inline injectate probes and kits compatible with BD or Edwards Lifesciences



Order Code	Description	Length/Size	Qty
<b>Cardiac Output Cables &amp; Probes</b>			
<b>Cables</b>			
 700-0027-00	TruLink Cardiac Output Adapter Cable	30 cm/12 in.	1 Each
 306655-003	Cardiac Output cable compatible with Baxter Edwards/Abbott injectate temperature probes and most models of thermodilution catheters		1 Each
 306655-004	Cardiac Output cable compatible with Becton Dickinson/Datex-Ohmeda and Argon Medical injectate temperature probes		1 Each
<b>Injectate Probes</b>			
 690-0023-00	In-line injectate temperature probe Compatible with Datex-Ohmeda SP4045		1 Each
 690-0024-00	Reference solution in-line injectate temperature probe. Compatible with Datex-Ohmeda SP5030		1 Each
 690-0033-00	In-line injectate temperature probe. Compatible with Datex-Ohmeda SP4042. Not compatible with Abbott or Baxter in-line injectate probes or delivery systems.		1 Each
 690-0043-01	In-line injectate sensor housing, single use, sterile. Compatible with Datex-Ohmeda SP5045		25/Box

## Cardiac Output Injectate Kits

### Closed Loop Injectate Kits

690-0031-00



CitiKit SP4572 CO in-line room temperature injectate kit, closed-loop delivery system, single use, sterile

5/Box

Indicates discontinued part, available while supplies last