

# SUPPLIES & ACCESSORIES


## ELECTROCARDIOGRAPHY

IEC Color Coding



### ECG Combiner Lead Wires & Cables

Hospitals are equipped with numerous monitors and other devices that use ECG lead wires, each with their own unique connectors which complicates patient transfers, undermines efficiency, and makes purchasing and maintenance more complex. Our ECG combiner lead wire system can standardize most equipment to one lead wire. The reusable or disposable lead wire sets are available in a variety of leads and lengths with either a snap or a multi-pin electrode interface. Validated for use with Spacelabs command modules 90496/91496 and Aria/Ultraview Telemetry systems.



Order Code	Description	Length	Qty
------------	-------------	--------	-----

#### ECG Combiner Lead Wire Sets - requires ECG Combiner Cable or Adapter



Multi-pin

#### Disposable, DIN Shielded, Multi-Pinch, IEC

700-0006-44	5 Lead	Multi-Pinch	127cm /50 in	25 Sets
700-0006-53	5 Lead (V-Lead)	Multi-Pinch	127cm /50 in	25 Sets



Multi-pin

#### Reusable, DIN Shielded, Multi-Pinch, IEC

700-0006-41	3 Lead	Multi-Pinch	74cm /29 in	1 Set
700-0006-43	5 Lead	Multi-Pinch	74cm /29 in	1 Set
700-0010-02	5 Lead (V-Lead)	Multi-Pinch	130cm/51 in	1 Set



Snap

#### Reusable, DIN Shielded, Snap, IEC

700-0007-41	3 Lead	Snap	74cm /29 in	1 Set
700-0007-44	5 Lead	Snap	74cm /29 in	1 Set

Indicates discontinued part, available while supplies last

#### ECG Combiner Cables

#### Spacelabs Connect Hardwired Cables for monitors using Command Module (90496/91496)



700-0008-81	3 Lead	IEC	305 cm/10 ft	1 Cable
-------------	--------	-----	--------------	---------



700-0008-83	5 Lead	IEC	305 cm/10 ft	1 Cable
-------------	--------	-----	--------------	---------



700-0008-84	10 Lead (5+5)	IEC	305 cm/10 ft	1 Cable
-------------	---------------	-----	--------------	---------

## ECG Lead Wires & Cables - Shielded

IEC color coded shielded ECG lead wires and cables for use with Spacelabs Command Modules 90496 & 91496 and mCare Monitor 91220.



Shielded Lead Wire



91496

Order Code	Description	Length	Qty
------------	-------------	--------	-----

### ECG Lead Wire Sets, Shielded, Reusable, IEC



Snap

#### Snap Sets

700-0007-01	3 Lead, shielded	61 cm/24 in	1 Set
700-0007-03	3 Lead, shielded	102 cm/40 in	1 Set
700-0007-09	5 Lead, shielded	61 cm/24 in	1 Set
700-0007-11	5 Lead, shielded	102 cm/40 in	1 Set
700-0007-17	10 Lead, limb/chest, shielded	61 cm/24 in	1 Set
700-0007-19	10 Lead, limb/chest, shielded	102 cm/40 in	1 Set



Pinch

#### Pinch Sets

700-0006-01	3 Lead, shielded	61 cm/24 in	1 Set
700-0006-03	3 Lead, shielded	102 cm/40 in	1 Set
700-0006-09	5 Lead, shielded	61 cm/24 in	1 Set
700-0006-11	5 Lead, shielded	102 cm/40 in	1 Set
700-0006-17	10 Lead, limb/chest, shielded	61 cm/24 in	1 Set
700-0006-19	10 Lead, limb/chest, shielded	102 cm/40 in	1 Set

### ECG Lead Wires, Individual Colors, Shielded, Reusable, IEC

#### ECG Cables, Shielded, IEC



#### With 17 PIN Adapter

700-0008-09	3 Lead, shielded	305 cm/10 ft	1 Cable
700-0008-07	5 Lead, shielded	305 cm/10 ft	1 Cable
700-0008-01	10 Lead, shielded	305 cm/10 ft	1 Cable

Indicates discontinued part, available while supplies last

#### ECG Accessories



#### Electrodes

690-0289-00	ECG electrode, transparent tape, adult, single use	Adult	30/Box
690-0262-00	ECG electrode, transparent tape, general/stress	Adult	30/Box



## ECG Lead Wires & Cables - Shielded




IEC color coded shielded ECG lead wires and cables for use with Spacelabs Command Modules 90496 & 91496 and mCare Monitor 91220.



Shielded Lead Wire



91496

Order Code	Description	Length	Qty
<b>Neonatal ECG Supplies &amp; Accessories</b>			
	<b>DIN Unshielded, Neonatal ECG Cable</b>		
700-0008-13	3 Lead		1 Cable
	<b>DIN Unshielded, Neonatal Electrodes</b>		
717-0080-00	Pre-wired neonatal electrodes, pre-gelled, silver/silver chloride, CE	3-Lead	100 Sets/Box
	<b>Electrode</b>		
690-0290-00	ECG electrode, transparent tape, pediatric/neonatal, single use	Ped/Neo	30/Box




Indicates discontinued part, available while supplies last

## DIN Unshielded Lead Wires & Telemetry

TruLink ECG Din unshielded reusable lead wires are validated for use with Spacelabs telemetry systems. They are available in 3 and 5-lead sets, or available to purchase individually in color-coded multi-packs. They can be ordered with a snap or pinch electrode interface. A grouper adapter is also available to help manage leads for a quick connection to the transmitter. Validated for use with models 90341, 90343, 90347, 91341, 91343, 91347 and 96281.



DIN Un-Shielded

Order Code	Description	Length/Size	Qty
<b>ECG Lead Wire Sets, DIN Unshielded, Reusable, IEC</b>			
	<b>Snap Sets</b>		
012-0328-10	5 Lead	61 cm/24 in	1 Set
	<b>Snap</b>		
	<b>Pinch Sets</b>		
012-0325-10	3 Lead	61 cm/24 in	1 Set
012-0327-10	5 Lead	61 cm/24 in	1 Set
	<b>Pinch</b>		
	<b>For use with DIN Unshielded Lead Wires</b>		
706-0190-00	Aria Telemetry Grouper (96281), IEC		5/Bag

### Telemetry Accessories



<b>Disposable Pouches and Electrodes</b>			
015-0500-02	Ultraview Telemetry pouch (90347)		100/Box
690-0289-00	ECG electrode, transparent tape, adult, single use	Adult	30/Box
690-0262-00	ECG transparent electrode, general/stress	Adult	30/Box









Indicates discontinued part, available while supplies last

## DIN Unshielded Lead Wires & Telemetry

TruLink ECG Din unshielded reusable lead wires are validated for use with Spacelabs telemetry systems. They are available in 3 and 5-lead sets, or available to purchase individually in color-coded multi-packs. They can be ordered with a snap or pinch electrode interface. A grouper adapter is also available to help manage leads for a quick connection to the transmitter. Validated for use with models 90341, 90343, 90347, 91341, 91343, 91347 and 96281.



DIN Un-Shielded

Order Code	Description	Length/Size	Qty
<b>SpO<sub>2</sub> Cable &amp; Sensors</b>			
<b>SpO<sub>2</sub> Cable</b>			
	700-0014-00 TruLink SpO2 Sensor Adapter Cable Compatible with SL multi-parameter telemetry transmitters (90343/91343/96281) and TruLink SpO2 sensors	60 cm/24 in	1 Cable
<b>Note: Only the below genuine Spacelabs and Nellcor branded SpO<sub>2</sub> sensors are compatible with Spacelabs Telemetry systems.</b>			
<b>Reusable SpO<sub>2</sub> Sensors</b>			
	015-0660-00 TruLink SpO <sub>2</sub> finger sensor	Adult >30 kg/66 lbs.	1 Each
	690-0003-01 Nellcor DS-100A, finger SpO2 sensor <i>Available for sale with monitor purchase</i>	Adult >30 kg/66 lbs.	1 Each
<b>Disposable SpO<sub>2</sub> Sensors</b>			
	015-0662-00 TruLink SpO <sub>2</sub> foam sensor	Adult > 30 kg/66 lbs.	24/Box
	015-0663-00 TruLink SpO <sub>2</sub> vinyl sensor	Adult >30 kg/66 lbs.	24/Box
	015-0664-00 TruLink SpO <sub>2</sub> foam sensor	Pediatric 10-50 kg/ 22-110 lbs.	24/Box
	015-0665-00 <b>TruLink SpO<sub>2</sub> vinyl sensor</b>	<b>Pediatric 10-50 kg/ 22-110 lbs.</b>	<b>24/Box</b>
	015-0667-00 TruLink SpO <sub>2</sub> cloth sensor	Neonatal <3 kg/7 lbs.	24/Box

Indicates discontinued part, available while supplies last