


SUPPLIES & ACCESSORIES









INVASIVE BLOOD PRESSURE



Invasive Pressure Supplies

Supplies are available through Spacelabs for 90496/91496 modules configured for invasive pressure. Customers can choose from Edwards Lifesciences or Abbott Laboratories integrated single or dual pressure cables that connect directly to the transducer(s). For BD models, Spacelabs offers single or dual adapters and interface cables to allow use of compatible transducers. Transducer kits are available in patient or pole mount options, as well as a neonatal configuration.



Order Code	Description		Length/Size	Qty
Integrated Adapter Cables (connects directly to transducer)				
Edward Lifesciences				
 700-0293-01	Single pressure integrated cable		396 cm/13 ft	1 Each
700-0292-01	Dual pressure integrated cable		396 cm/13 ft	1 Each
Abbott Laboratories/Hospira/ICU Medical				
 700-0295-01	Single pressure integrated cable		396 cm/13 ft	1 Each
700-0294-01	Dual pressure integrated cable		396 cm/13 ft	1 Each
Interface Cable & Adapters				
Interface Cable (Datex-Ohmeda TC-VTK or Becton Dickinson model 684085 pressure transducers)				
 690-0021-00	Requires adapter cable 700-0091-00 or 700-0028-00		366 cm/12 ft	1 Each
Interface Adapter- Transducer interface cable required from transducer company				
 700-0791-00	Single pressure cable		122 cm/4 ft	1 Each
700-0028-00	Dual pressure cable		122 cm/4 ft	1 Each
Transducer Kits & Accessories				
Pressure Transducer Kits, Single use				
 690-0008-00	Transducer Kit, Patient Mount, DTX/Plus		152.5 cm/60 in	5/Box
 690-0009-00	Transducer Kit, Pole Mount, DTX/Plus		152.5 cm/60 in	5/Box

Indicates discontinued part, available while supplies last