

## TruLink® Nylon Reusable

TruLink reusable nylon cuffs provide an economical and durable solution for your NIBP measurement needs. TruLink cuffs feature an anti-microbial coating, fully cleanable nylon construction, and are durable for multi-patient extended use. Available in 8 sizes, TruLink cuffs provide color-coded sizing and markings to assist you in selecting the right size for you patient. TruLink cuffs are naturally latex and PVC-free.



Range/  
Length Qty

Order Code Description Size

### TruLink Nylon Reusable BP Cuffs

#### Single Tube

*For use with Command Modules 90496/91496 and Ultraview DM3 (91330) and DM4 (91331) Monitors*

715-1111-10	Single tube, HP	Infant	8 - 14 cm	1 Each
715-1121-10	Single tube, HP	Child	13 - 20 cm	1 Each
715-1131-10	Single tube, HP	Small Adult	18 - 26 cm	1 Each
715-1141-10	Single tube, HP	Adult	26 - 35 cm	1 Each
715-1151-10	Single tube, HP	Adult Long	29 - 38 cm	1 Each
715-1161-10	Single tube, HP	Large Adult	32 - 42 cm	1 Each
715-1171-10	Single tube, HP	Large Adult Long	35 - 44 cm	1 Each
715-1181-10	Single tube, HP	Thigh Cuff	42 - 50 cm	1 Each



Bayonet



#### Single Tube Hose

714-0018-00	Single tube, quick connect	Adult	275 cm/9 ft	1 Each
714-0060-00	Single tube, quick connect	Adult	366 cm/12 ft	1 Each

## TruLink® Soft Disposable

TruLink fabric cuffs provide an economical and patient friendly solution for your NIBP measurement needs. TruLink cuffs feature a super-soft Sontara fabric which is comfortable, helps reduce perspiration under the cuff, and is durable for long patient stays. The single patient - multiple use application offers a cost effective solution for infection control. Available in 8 sizes, TruLink cuffs provide a wide range of options for correct readings.



Order Code	Description	Size	Range/ Length	Qty
------------	-------------	------	------------------	-----

### TruLink Soft Disposable BP Cuffs

#### Single Tube

*For use with Command Modules 90496/91496 and Ultraview DM3 (91330) and DM4 (91331) Monitors*

715-1113-10	Single tube, HP	Infant	8 - 14 cm	5/Bag
715-1123-10	Single tube, HP	Child	13 - 20 cm	5/Bag
715-1133-10	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
715-1143-10	Single tube, HP	Adult	26 - 35 cm	5/Bag
715-1153-10	Single tube, HP	Adult Long	29 - 38 cm	5/Bag
715-1163-10	Single tube, HP	Large Adult	32 - 42 cm	5/Bag
715-1173-10	Single tube, HP	Large Adult Long	35 - 44 cm	5/Bag
715-1183-10	Single tube, HP	Thigh Cuff	42 - 50 cm	5/Bag



Bayonet

#### Single Tube Hose

714-0018-00	Single tube, quick connect	Adult	275 cm/9 ft	1 Each
714-0060-00	Single tube, quick connect	Adult	366 cm/12 ft	1 Each



## TruLink® Vinyl Disposable

TruLink vinyl cuffs provide an economical and patient friendly solution for your NIBP measurement needs. TruLink cuffs feature a soft vinyl construction which is comfortable, cleanable, and is durable for long patient stays. The single patient - multiple use application offers a cost effective solution for infection control. Available in 8 sizes, TruLink cuffs provide a wide range of options for correct readings.



Order Code	Description	Size	Range/ Length	Qty
------------	-------------	------	------------------	-----

### TruLink Vinyl Disposable BP Cuffs

#### Single Tube

*For use with Command Modules 90496/91496 and Ultraview DM3 (91330) and DM4 (91331) Monitors*

715-1112-10	Single tube, HP	Infant	8 - 14 cm	5/Bag
715-1122-10	Single tube, HP	Child	13 - 20 cm	5/Bag
715-1132-10	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
715-1142-10	Single tube, HP	Adult	26 - 35 cm	5/Bag
715-1152-10	Single tube, HP	Adult Long	29 - 38 cm	5/Bag
715-1162-10	Single tube, HP	Large Adult	32 - 42 cm	5/Bag
715-1172-10	Single tube, HP	Large Adult Long	35 - 44 cm	5/Bag
715-1182-10	Single tube, HP	Thigh Cuff	42 - 50 cm	5/Bag



Bayonet

#### Single Tube Hose

714-0018-00	Single tube, quick connect	Adult	275 cm/9 ft	1 Each
714-0060-00	Single tube, quick connect	Adult	366 cm/12 ft	1 Each



## Blood Pressure Cuff Accessories

Spacelabs non-invasive blood pressure cuff accessories are compatible with all Spacelabs and TruLink branded non-invasive blood pressure cuffs. Choose from a variety of hoses, adapters and conversion accessories which offer ease-of-use and flexibility.


Order Code	Description	Size	Length	Qty
<b>Non-Invasive Blood Pressure Hose</b>				
<b>Single Tube Hose</b>				
<i>For use with Command Modules 90496/91496 and Ultraview DM3 (91330) and DM4 (91331) Monitors</i>				
714-0018-00	90496/91496/91330	Adult	275 cm/9 ft	1 Each
714-0060-00	90496/91496/91330	Adult	366 cm/12 ft	1 Each
714-0019-01	90496/91496	Neonatal	275 cm/9 ft	1 Each
<b>Double Tube to Single Tube Hose</b>				
<i>For use with Command Modules 90496/91496</i>				
714-0018-12	90496/91496	Adult	275 cm/9 ft	1 Each
714-0019-12	90496/91496	Neonatal	275 cm/9 ft	1 Each



Indicates discontinued part, available while supplies last

## Hose Connectors & Adapters





<b>Slip Lock Fitting</b>				
<i>Fitting for Spacelabs standard NIBP tubing. Adapts hoses for use with quick connect fittings.</i>				
714-0038-00	3mm ID Adult NIBP hoses			1 Each
714-0047-00	4mm ID Adult NIBP hoses			1 Each
714-0046-00	5mm ID Adult NIBP hoses			1 Each
714-0080-00	Dual male screw to single female bayonet Y adapter	Attaches HP bayonet cuff to dual screw hose		1 Each
714-0081-00	Dual female screw to single male bayonet Y adapter	Attaches dual screw cuff to Spacelabs hose (714-0018-00)		1 Each
714-0082-00	MQML to single male bayonet Y adapter	Attaches MQML cuff to Spacelabs hose (714-0018-00)		1 Each
714-0083-00	MQML to single female bayonet Y adapter	Attaches HP bayonet cuff to MQML hose		1 Each
714-0084-00	Dual MQ to single female bayonet Y adapter	Attaches Spacelabs cuff to dual MQ hose		1 Each



## UltraCheck® Nylon Reusable

UltraCheck reusable nylon cuffs provide an economical and durable solution for your NIBP measurement needs. UltraCheck cuffs feature an anti-microbial coating, fully cleanable nylon construction, and are durable for multi-patient extended use. Available in 8 sizes, UltraCheck cuffs provide color-coded sizing and markings to assist you in selecting the right size for you patient. UltraCheck cuffs are naturally latex and PVC-free.






	Order Code	Description	Size	Length	Qty
<b>UltraCheck Nylon Reusable BP Cuffs</b>					
<b>Single Tube</b>					
 Plastic Bayonet	US0814HP-05	Single tube, HP	Infant	8 - 14 cm	5/Bag
	US1320HP-05	Single tube, HP	Child	13 - 20 cm	5/Bag
	US1826HP-05	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
	US2635HP-05	Single tube, HP	Adult	26 - 35 cm	5/Bag
	US2938HP-05	Single tube, HP	Adult Long	29 - 38 cm	5/Bag
	US3242HP-05	Single tube, HP	Large Adult	32 - 42 cm	5/Bag
	US3544HP-05	Single tube, HP	Large Adult Long	35 - 44 cm	5/Bag
	US4250HP-05	Single tube, HP	Thigh	42 - 50 cm	5/Bag
 Metal Bayonet	US0814MHP-05	Single tube, MHP	Infant	8 - 14 cm	5/Bag
	US1320MHP-05	Single tube, MHP	Child	13 - 20 cm	5/Bag
	US1826MHP-05	Single tube, MHP	Small Adult	18 - 26 cm	5/Bag
	US2635MHP-05	Single tube, MHP	Adult	26 - 35 cm	5/Bag
	US2938MHP-05	Single tube, MHP	Adult Long	29 - 38 cm	5/Bag
	US3242MHP-05	Single tube, MHP	Large Adult	32 - 42 cm	5/Bag
	US3544MHP-05	Single tube, MHP	Large Adult Long	35 - 44 cm	5/Bag
	US4250MHP-05	Single tube, MHP	Thigh	42 - 50 cm	5/Bag
 Threaded Screw	US1320DM-05	Single tube, DM	Child	13 - 20 cm	5/Bag
	US1826DM-05	Single tube, DM	Small Adult	18 - 26 cm	5/Bag
	US2635DM-05	Single tube, DM	Adult	26 - 35 cm	5/Bag
	US2938DM-05	Single tube, DM	Adult Long	29 - 38 cm	5/Bag
	US3544DM-05	Single tube, DM	Large Adult	32 - 42 cm	5/Bag
	US3242DM-05	Single tube, DM	Large Adult	32 - 42 cm	5/Bag
	US4250DM-05	Single tube, DM	Thigh	42 - 50 cm	5/Bag
	 Female Subminiature Twist Lock	US1320MQ-05	Single tube, MQ	Child	13 - 20 cm
US1826MQ-05		Single tube, MQ	Small Adult	18 - 26 cm	5/Bag
US2635MQ-05		Single tube, MQ	Adult	26 - 35 cm	5/Bag
US2938MQ-05		Single tube, MQ	Adult Long	29 - 38 cm	5/Bag
US3242MQ-05		Single tube, MQ	Large Adult	32 - 42 cm	5/Bag
US3544MQ-05		Single tube, MQ	Large Adult Long	35 - 44 cm	5/Bag
US4250MQ-05		Single tube, MQ	Thigh	42 - 50 cm	5/Bag

Indicates discontinued part, available while supplies last

## UltraCheck® Nylon Reusable

UltraCheck reusable nylon cuffs provide an economical and durable solution for your NIBP measurement needs. UltraCheck cuffs feature an anti-microbial coating, fully cleanable nylon construction, and are durable for multi-patient extended use. Available in 8 sizes, UltraCheck cuffs provide color-coded sizing and markings to assist you in selecting the right size for you patient. UltraCheck cuffs are naturally latex and PVC-free.







	Order Code	Description	Size	Length	Qty
<b>UltraCheck Nylon Reusable BP Cuffs</b>					
<b>Single Tube</b>					
 <b>Male Subminiature Twist Lock</b>	US0814MQML-05	Single tube, MQML	Infant	8 - 14 cm	5/Bag
	US1320MQML-05	Single tube, MQML	Child	13 - 20 cm	5/Bag
	US1826MQML-05	Single tube, MQML	Small Adult	18 - 26 cm	5/Bag
	US2635MQML-05	Single tube, MQML	Adult	26 - 35 cm	5/Bag
	US2938MQML-05	Single tube, MQML	Adult Long	29 - 38 cm	5/Bag
	US3242MQML-05	Single tube, MQML	Large Adult	32 - 42 cm	5/Bag
	US3544MQML-05	Single tube, MQML	Large Adult Long	35 - 44 cm	5/Bag
 <b>No Fitting</b>	US0814-05	Single tube, no fitting	Infant	8 - 14 cm	5/Bag
	US1320-05	Single tube, no fitting	Child	13 - 20 cm	5/Bag
	US1826-05	Single tube, no fitting	Small Adult	18 - 26 cm	5/Bag
	US2635-05	Single tube, no fitting	Adult	26 - 35 cm	5/Bag
	US2938-05	Single tube, no fitting	Adult Long	29 - 38 cm	5/Bag
	US3242-05	Single tube, no fitting	Large Adult	32 - 42 cm	5/Bag
	US3544-05	Single tube, no fitting	Large Adult Long	35 - 44 cm	5/Bag
 <b>Manual Bulb</b>	US0814MB	Single tube, MB w/gauge	Infant	8 - 14 cm	1 Each
	US1826MB	Single tube, MB w/gauge	Small Adult	18 - 26 cm	1 Each
	US2635MB	Single tube, MB w/gauge	Adult	26 - 35 cm	1 Each
	US2938MB	Single tube, MB w/gauge	Adult Long	29 - 38 cm	1 Each
	US3242MB	Single tube, MB w/gauge	Large Adult	32 - 42 cm	1 Each
	US3544MB	Single tube, MB w/gauge	Large Adult Long	35 - 44 cm	1 Each
	US4250MB	Single tube, MB w/gauge	Thigh	42 - 50 cm	1 Each
Indicates discontinued part, available while supplies last					

## UltraCheck® Nylon Reusable

UltraCheck reusable nylon cuffs provide an economical and durable solution for your NIBP measurement needs. UltraCheck cuffs feature an anti-microbial coating, fully cleanable nylon construction, and are durable for multi-patient extended use. Available in 8 sizes, UltraCheck cuffs provide color-coded sizing and markings to assist you in selecting the right size for you patient. UltraCheck cuffs are naturally latex and PVC-free.



	Order Code	Description	Size	Length	Qty
<b>UltraCheck Nylon Reusable BP Cuffs</b>					
<b>Double Tube</b>					
 <b>Threaded Screw</b>	UD0814DM-05	Double tube, DM	Infant	8 - 14 cm	5/Bag
	UD1320DM-05	Double tube, DM	Child	13 - 20 cm	5/Bag
	UD1826DM-05	Double tube, DM	Small Adult	18 - 26 cm	5/Bag
	UD2635DM-05	Double tube, DM	Adult	26 - 35 cm	5/Bag
	UD2938DM-05	Double tube, DM	Adult Long	29 - 38 cm	5/Bag
	UD3242DM-05	Double tube, DM	Large Adult	32 - 42 cm	5/Bag
	UD3544DM-05	Double tube, DM	Large Adult Long	35 - 44 cm	5/Bag
	UD4250DM-05	Double tube, DM	Thigh	42 - 50 cm	5/Bag
 <b>Female Subminiature Twist Lock</b>	UD0814MQ-05	Double tube, MQ	Infant	8 - 14 cm	5/Bag
	UD1320MQ-05	Double tube, MQ	Child	13 - 20 cm	5/Bag
	UD1826MQ-05	Double tube, MQ	Small Adult	18 - 26 cm	5/Bag
	UD2635MQ-05	Double tube, MQ	Adult	26 - 35 cm	5/Bag
	UD2938MQ-05	Double tube, MQ	Adult Long	29 - 38 cm	5/Bag
	UD3242MQ-05	Double tube, MQ	Large Adult	32 - 42 cm	5/Bag
	UD3544MQ-05	Double tube, MQ	Large Adult Long	35 - 44 cm	5/Bag
	UD4250MQ-05	Double tube, MQ	Thigh	42 - 50 cm	5/Bag
 <b>Male &amp; Female Mated</b>	UD0814MQMF-05	Double tube, MQMF	Infant	8 - 14 cm	5/Bag
	UD1320MQMF-05	Double tube, MQMF	Child	13 - 20 cm	5/Bag
	UD1826MQMF-05	Double tube, MQMF	Small Adult	18 - 26 cm	5/Bag
	UD2635MQMF-05	Double tube, MQMF	Adult	26 - 35 cm	5/Bag
	UD2938MQMF-05	Double tube, MQMF	Adult Long	29 - 38 cm	5/Bag
	UD3242MQMF-05	Double tube, MQMF	Large Adult	32 - 42 cm	5/Bag
	UD4250MQMF-05	Double tube, MQMF	Thigh	42 - 50 cm	5/Bag
	 <b>No Fitting</b>	UD0814-05	Double tube, no fitting	Infant	8 - 14 cm
UD1320-05		Double tube, no fitting	Child	13 - 20 cm	5/Bag
UD1826-05		Double tube, no fitting	Small Adult	18 - 26 cm	5/Bag
UD2635-05		Double tube, no fitting	Adult	26 - 35 cm	5/Bag
UD2938-05		Double tube, no fitting	Adult Long	29 - 38 cm	5/Bag
UD3242-05		Double tube, no fitting	Large Adult	32 - 42 cm	5/Bag
UD3544-05		Double tube, no fitting	Large Adult Long	35 - 44 cm	5/Bag
UD4250-05		Double tube, no fitting	Thigh	42 - 50 cm	5/Bag

## UltraCheck® Nylon Reusable

UltraCheck reusable nylon cuffs provide an economical and durable solution for your NIBP measurement needs. UltraCheck cuffs feature an anti-microbial coating, fully cleanable nylon construction, and are durable for multi-patient extended use. Available in 8 sizes, UltraCheck cuffs provide color-coded sizing and markings to assist you in selecting the right size for you patient. UltraCheck cuffs are naturally latex and PVC-free.



Order Code	Description	Size	Length	Qty
------------	-------------	------	--------	-----

### UltraCheck Nylon Reusable BP Cuffs

#### Double Tube



UD0814MB-05	Double tube, MB	Infant	8 - 14 cm	5/Bag
UD1826MB-05	Double tube, MB	Small Adult	18 - 26 cm	5/Bag
UD2635MB-05	Double tube, MB	Adult	26 - 35 cm	5/Bag
UD2938MB-05	Double tube, MB	Adult Long	29 - 38 cm	5/Bag
UD3242MB-05	Double tube, MB	Large Adult	32 - 42 cm	5/Bag
UD3544MB-05	Double tube, MB	Large Adult Long	35 - 44 cm	5/Bag
UD4250MB-05	Double tube, MB	Thigh	42 - 50 cm	5/Bag

Indicates discontinued part, available while supplies last



## SoftCheck® Fabric Disposable

SoftCheck fabric cuffs provide an economical and patient friendly solution for your NIBP measurement needs. SoftCheck cuffs feature a super-soft Sontara fabric which is comfortable, helps reduce perspiration under the cuff, and is durable for long patient stays. The single patient - multiple use application offers a cost effective solution for infection control. Available in 8 sizes, SoftCheck cuffs provide a wide range of options for correct readings.



Range/  
Length

Qty

Order Code

Description

Size

### SoftCheck Fabric Disposable BP Cuffs

#### Single Tube



Plastic Bayonet

ST0814HP-05	Single tube, HP	Infant	8 - 14 cm	5/Bag
<b>ST0814HP-20</b>	<b>Single tube, HP</b>	<b>Infant</b>	<b>8 - 14 cm</b>	<b>20/Box</b>
ST1320HP-05	Single tube, HP	Child	13 - 20 cm	5/Bag
ST1320HP-20	Single tube, HP	Child	13 - 20 cm	20/Box
ST1826HP-05	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
ST1826HP-20	Single tube, HP	Small Adult	18 - 26 cm	20/Box
ST2635HP-05	Single tube, HP	Adult	26 - 35 cm	5/Bag
ST2635HP-20	Single tube, HP	Adult	26 - 35 cm	20/Box
ST2938HP-05	Single tube, HP	Adult Long	29 - 38 cm	5/Bag
ST2938HP-20	Single tube, HP	Adult Long	29 - 38 cm	20/Box
ST3242HP-05	Single tube, HP	Large Adult	32 - 42 cm	5/Bag
ST3242HP-20	Single tube, HP	Large Adult	32 - 42 cm	20/Box
ST3544HP-05	Single tube, HP	Large Adult Long	35 - 44 cm	5/Bag
ST3544HP-20	Single tube, HP	Large Adult Long	35 - 44 cm	20/Box
ST4250HP-05	Single tube, HP	Thigh	42 - 50 cm	5/Bag
<b>ST4250HP-20</b>	<b>Single tube, HP</b>	<b>Thigh</b>	<b>42 - 50 cm</b>	<b>20/Box</b>



Female  
Subminiature

<b>ST0814MQ-05</b>	<b>Single tube, MQ</b>	<b>Infant</b>	<b>8 - 14 cm</b>	<b>5/Bag</b>
ST1320MQ-05	Single tube, MQ	Child	13 - 20 cm	5/Bag
ST1826MQ-05	Single tube, MQ	Small Adult	18 - 26 cm	5/Bag
ST2635MQ-05	Single tube, MQ	Adult	26 - 35 cm	5/Bag
ST2938MQ-05	Single tube, MQ	Adult Long	29 - 38 cm	5/Bag
ST3242MQ-05	Single tube, MQ	Large Adult	32 - 42 cm	5/Bag



Threaded Screw

<b>ST1826DM-05</b>	<b>Single tube, DM</b>	<b>Small Adult</b>	<b>18 - 26 cm</b>	<b>5/Bag</b>
<b>ST2635DM-05</b>	<b>Single tube, DM</b>	<b>Adult</b>	<b>26 - 35 cm</b>	<b>5/Bag</b>
<b>ST2938DM-05</b>	<b>Single tube, DM</b>	<b>Adult Long</b>	<b>29 - 38 cm</b>	<b>5/Bag</b>
<b>ST3242DM-05</b>	<b>Single tube, DM</b>	<b>Large Adult</b>	<b>32 - 42 cm</b>	<b>5/Bag</b>

Indicates discontinued part, available while supplies last

## SoftCheck® Fabric Disposable

SoftCheck fabric cuffs provide an economical and patient friendly solution for your NIBP measurement needs. SoftCheck cuffs feature a super-soft Sontara fabric which is comfortable, helps reduce perspiration under the cuff, and is durable for long patient stays. The single patient - multiple use application offers a cost effective solution for infection control. Available in 8 sizes, SoftCheck cuffs provide a wide range of options for correct readings.



Range/  
Length

Qty

Order Code

Description

Size

### SoftCheck Fabric Single Patient Multiple Use BP Cuffs



No Fitting

#### Single Tube

ST0814-05	Single tube, no fitting	Infant	8 - 14 cm	5/Bag
ST1320-05	Single tube, no fitting	Child	13 - 20 cm	5/Bag
ST2635-05	Single tube, no fitting	Adult	26 - 35 cm	5/Bag
ST3544-05	Single tube, no fitting	Large Adult Long	35 - 44 cm	5/Bag
ST4250-05	Single tube, no fitting	Thigh	42 - 50 cm	5/Bag



Threaded Screw

#### Double Tube

DT2635DM-05	Double tube, DM	Adult	26 - 35 cm	5/Bag
DT2635DM-20	Double tube, DM	Adult	26 - 35 cm	20/Box



Twist Lock

DT2635MQ-20	Double tube, MQ	Adult	26 - 35 cm	20/Box
-------------	-----------------	-------	------------	--------



No fitting

DT2635MB-05	Double tube, MB	Adult	26 - 35 cm	5/Box
-------------	-----------------	-------	------------	-------

DT2938-05	Double tube, no fitting	Adult Long	29 - 38 cm	5/Bag
-----------	-------------------------	------------	------------	-------

Indicates discontinued part, available while supplies last

## SoftCheck® Vinyl Disposable

SoftCheck vinyl cuffs provide an economical and patient friendly solution for your NIBP measurement needs. SoftCheck cuffs feature a soft vinyl construction which is comfortable, cleanable, and is durable for long patient stays. The single patient - multiple use application offers a cost effective solution for infection control. Available in 8 sizes, SoftCheck cuffs provide a wide range of options for correct readings.



Order Code	Description	Size	Range/ Length	Qty
------------	-------------	------	------------------	-----

### SoftCheck Vinyl Disposables BP Cuffs

#### Single Tube



Plastic Bayonet

VST0814HP-05	Single tube, HP	Infant	8 - 14 cm	5/Bag
VST1320HP-05	Single tube, HP	Child	13 - 20 cm	5/Bag
VST1826HP-05	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
VST2635HP-05	Single tube, HP	Adult	26 - 35 cm	5/Bag
VST2635HP-20	Single tube, HP	Adult	26 - 35 cm	20/Box
VST2938HP-05	Single tube, HP	Adult Long	29 - 38 cm	5/Bag
VST3242HP-05	Single tube, HP	Large Adult	32 - 42 cm	5/Bag
VST3544HP-05	Single tube, HP	Large Adult Long	35 - 44 cm	5/Bag
VST4250HP-05	Single tube, HP	Thigh	42 - 50 cm	5/Bag



Subminiature  
Twist Lock

VST0814MQ-05	Single tube, MQ	Infant	8 - 14 cm	5/Bag
VST1320MQ-05	Single tube, MQ	Child	13 - 20 cm	5/Bag
VST1826MQ-05	Single tube, MQ	Small Adult	18 - 26 cm	5/Bag
VST2635MQ-20	Single tube, MQ	Adult	26 - 35 cm	20/Box
VST3242MQ-05	Single tube, MQ	Large Adult	32 - 42 cm	5/Bag
VST4250MQ-05	Single tube, MQ	Thigh	42 - 50 cm	5/Bag



Male  
Subminiature  
Twist Lock

VST0814MQML5-20	Single tube, MQ	Infant	8 - 14 cm	5/Bag
VST1320MQML5-20	Single tube, MQ	Child	13 - 20 cm	5/Bag
VST1826MQML5-20	Single tube, MQ	Small Adult	18 - 26 cm	5/Bag
VST2635MQML5-20	Single tube, MQ	Adult	26 - 35 cm	5/Bag
VST2938MQML5-20	Single tube, MQ	Adult Long	29 - 38 cm	20/Box

Indicates discontinued part, available while supplies last

**UltraCheck® Curve Nylon Reusable**  
**SoftCheck® Curve Fabric Disposable**

Your patients come in all shapes and sizes. That's why we offer the Curve cuff for extra-large adults. Its conical-shaped design provides a better fit for the contours of your patient's arm. This can increase patient comfort while providing a more accurate reading.



Order Code	Description	Size	Range/ Length	Qty
------------	-------------	------	------------------	-----

**UltraCheck Curve Nylon Reusable BP Cuff**

Single Tube				
BRUS3854	Single tube, no fitting	Extra Large Adult	38 - 54 cm	1 Each
BRUS3854HP	Single tube, plastic HP	Extra Large Adult	38 - 54 cm	1 Each
BRUS3854MHP	Single tube, metal HP	Extra Large Adult	38 - 54 cm	1 Each
BRUS3854MB	Single tube, manual bulb	Extra Large Adult	38 - 54 cm	1 Each
Double Tube				
BRUD3854	Double tube, no fitting	Extra Large Adult	38 - 54 cm	1 Each
BRUD3854DM	Double tube, DM	Extra Large Adult	38 - 54 cm	1 Each
BRUD3854MQ	Double tube, MQ	Extra Large Adult	38 - 54 cm	1 Each
BRUD3854MQMF	Double tube, MQMF	Extra Large Adult	38 - 54 cm	1 Each
BRUD3854MQML	Double tube, MQML	Extra Large Adult	38 - 54 cm	1 Each
BRUD3854MB	Double tube, manual bulb	Extra Large Adult	38 - 54 cm	1 Each

**SoftCheck Curve Fabric Disposable BP Cuff**

Single Tube				
BRST3854-05	Single tube, no fitting	Extra Large Adult	38 - 54 cm	5/Box
BRST3854HP-05	Single tube, plastic HP	Extra Large Adult	38 - 54 cm	5/Box
Double Tube				
BRDT3854DM-05	Double tube, DM	Extra Large Adult	38 - 54 cm	5/Box

Indicates discontinued part, available while supplies last



## SoftCheck® Fabric Neonatal Disposable

SoftCheck fabric cuffs provide an economical and patient friendly solution for your NIBP measurement needs. SoftCheck cuffs feature a super-soft Sontara fabric which is comfortable, helps reduce perspiration under the cuff, and is durable for long patient stays. The single patient - multiple use application offers a cost effective solution for infection control. Available in 5 different sizes for accurate readings.



Order Code	Description	Size	Range/ Length	Qty
------------	-------------	------	------------------	-----

### SoftCheck Fabric Neonatal Disposable BP Cuffs



Single Tube					
NN1ST-M316-10	Single tube, M316	Size 1	3-6 cm	10/Box	
NN2ST-M316-10	Single tube, M316	Size 2	4-8 cm	10/Box	
NN3ST-M316-10	Single tube, M316	Size 3	6-11 cm	10/Box	
NN4ST-M316-10	Single tube, M316	Size 4	7-13 cm	10/Box	
NN5ST-M316-10	Single tube, M316	Size 5	8-15 cm	10/Box	
NN1ST-10	Single tube, ML	Size 1	3-6 cm	10/Box	
NN2ST-10	Single tube, ML	Size 2	4-8 cm	10/Box	
NN3ST-10	Single tube, ML	Size 3	6-11 cm	10/Box	
NN4ST-10	Single tube, ML	Size 4	7-13 cm	10/Box	
NN5ST-10	Single tube, ML	Size 5	8-15 cm	10/Box	
NN1ST-HP-10	Single tube, HP	Size 1	3-6 cm	10/Box	
NN2ST-HP-10	Single tube, HP	Size 2	4-8 cm	10/Box	
NN5ST-HP-10	Single tube, HP	Size 5	8-15 cm	10/Box	



Indicates discontinued part, available while supplies last



Double Tube					
NN1DT-10	Double tube, ML	Size 1	3-6 cm	10/Box	
NN2DT-10	Double tube, ML	Size 2	4-8 cm	10/Box	
NN3DT-10	Double tube, ML	Size 3	6-11 cm	10/Box	
NN4DT-10	Double tube, ML	Size 4	7-13 cm	10/Box	
NN5DT-10	Double tube, ML	Size 5	8-15 cm	10/Box	






## SoftCheck® Vinyl Neonatal Disposable

SoftCheck vinyl cuffs provide an economical and patient friendly solution for your NIBP measurement needs. SoftCheck cuffs feature a soft vinyl construction which is comfortable, cleanable, and is durable for long patient stays. The single patient - multiple use application offers a cost effective solution for infection control. Available in 5 different sizes for accurate readings.



Order Code	Description	Size	Range/ Length	Qty
------------	-------------	------	------------------	-----

### SoftCheck Vinyl Neonatal Disposable BP Cuffs

Single Tube					
	VNN1ST-M316-10	Single tube, M316	Size 1	3-6 cm	10/Box
	VNN2ST-M316-10	Single tube, M316	Size 2	4-8 cm	10/Box
	VNN3ST-M316-10	Single tube, M316	Size 3	6-11 cm	10/Box
	VNN4ST-M316-10	Single tube, M316	Size 4	7-13 cm	10/Box
	VNN5ST-M316-10	Single tube, M316	Size 5	8-15 cm	10/Box
	VNN1ST-SQ-10	Single tube, SQ	Size 1	3-6 cm	10/Box
	VNN2ST-SQ-10	Single tube, SQ	Size 2	4-8 cm	10/Box
	VNN3ST-SQ-10	Single tube, SQ	Size 3	6-11 cm	10/Box
	VNN4ST-SQ-10	Single tube, SQ	Size 4	7-13 cm	10/Box
	VNN5ST-SQ-10	Single tube, SQ	Size 5	8-15 cm	10/Box
	VNN1ST-10	Single tube, ML	Size 1	3-6 cm	10/Box
	VNN2ST-10	Single tube, ML	Size 2	4-8 cm	10/Box
	VNN3ST-10	Single tube, ML	Size 3	6-11 cm	10/Box
	VNN4ST-10	Single tube, ML	Size 4	7-13 cm	10/Box
	VNN5ST-10	Single tube, ML	Size 5	8-15 cm	10/Box
	VNN1ST-HP-10	Single tube, HP	Size 1	3-6 cm	10/Box
	VNN5ST-HP-10	Single tube, HP	Size 5	8-15 cm	10/Box
	VNN1ST-MQ-10	Single tube, MQ	Size 1	3-6 cm	10/Box
	VNN2ST-MQ-10	Single tube, MQ	Size 2	4-8 cm	10/Box
	VNN3ST-MQ-10	Single tube, MQ	Size 3	6-11 cm	10/Box
	VNN4ST-MQ-10	Single tube, MQ	Size 4	7-13 cm	10/Box
	VNN5ST-MQ-10	Single tube, MQ	Size 5	8-15 cm	10/Box

Indicates discontinued part, available while supplies last

## SoftCheck® Vinyl Neonatal Disposable

SoftCheck vinyl cuffs provide an economical and patient friendly solution for your NIBP measurement needs. SoftCheck cuffs feature a soft vinyl construction which is comfortable, cleanable, and is durable for long patient stays. The single patient - multiple use application offers a cost effective solution for infection control. Available in 5 different sizes for accurate readings.



Order Code	Description	Size	Range/ Length	Qty
<b>Double Tube</b>				
VNN1DT-10	Double tube, ML	Size 1	3-6 cm	10/Box
VNN2DT-10	Double tube, ML	Size 2	4-8 cm	10/Box
VNN3DT-10	Double tube, ML	Size 3	6-11 cm	10/Box
VNN4DT-10	Double tube, ML	Size 4	7-13 cm	10/Box
VNN5DT-10	Double tube, ML	Size 5	8-15 cm	10/Box



Male Luer