

SUPPLIES & ACCESSORIES

PRESSURE INFUSORS



Unifusor® Pressure Infusors

Unifusor Pressure Infusors are available in 3 Models: Unifusor Disposable, Unifusor II & Unifusor Plus Reusable. All models are available in 500, 100 and 3,000 mL sizes. The Unifusors are 100% tested for leakage during the manufacturing process. Unifusor pressure Infusors include an IV pole loop, a fluid bag hood and a reinforced translucent sleeve which enables easy viewing of fluid type and level. This unit is designed for rapid IV infusion, invasive pressure monitoring and auto-transfusion procedures.



Order Code	Description	Size	Qty
------------	-------------	------	-----

Unifusor Disposable Pressure Infusor



Disposable Infusor with aneroid gauge - sold as one piece/complete unit

803FGA-24	Unifusor with thumbwheel valve & aneroid gauge	500mL	24/Box
804FGA-24	Unifusor with thumbwheel valve & aneroid gauge	1000mL	24/Box
805FGA-10	Unifusor with thumbwheel valve & aneroid gauge	3000mL	10/Box



803SGA-24	Unifusor with stopcock valve & aneroid gauge	500mL	24/Box
804SGA-24	Unifusor with stopcock valve & aneroid gauge	1000mL	24/Box
805SGA-10	Unifusor with stopcock valve & aneroid gauge	3000mL	10/Box



Disposable Infusor with piston gauge - sold as one piece/complete unit

903TGA-24	Unifusor with thumbwheel valve & piston gauge	500mL	24/Box
904TGA-24	Unifusor with thumbwheel valve & piston gauge	1000mL	24/Box
905TGA-10	Unifusor with thumbwheel valve & piston gauge	3000mL	10/Box



903SGA-24	Unifusor with stopcock valve & piston gauge	500mL	24/Box
904SGA-24	Unifusor with stopcock valve & piston gauge	1000mL	24/Box
905SGA-10	Unifusor with stopcock valve & piston gauge	3000mL	10/Box

Unifusor® Pressure Infusors

Unifusor Pressure Infusors are available in 3 Models: Unifusor Disposable, Unifusor II & Unifusor Plus Reusable. All models are available in 500, 100 and 3,000 mL sizes. The Unifusors are 100% tested for leakage during the manufacturing process. Unifusor pressure Infusors include an IV pole loop, a fluid bag hood and a reinforced translucent sleeve which enables easy viewing of fluid type and level. This unit is designed for rapid IV infusion, invasive pressure monitoring and auto-transfusion procedures.



Order Code	Description	Size	Qty
------------	-------------	------	-----

Unifusor II Disposable Pressure Infusor with Reusable Inflation System



Disposable Infusor - white disposable bag with reusable inflation system

2803FGA-24	Unifusor II with aneroid gauge & thumbwheel valve	500mL	24/Box
2804FGA-24	Unifusor II with aneroid gauge & thumbwheel valve	1000mL	24/Box
2805FGA-10	Unifusor II with aneroid gauge & thumbwheel valve	3000mL	10/Box



2803SGA-24	Unifusor II with aneroid gauge & stopcock valve	500mL	24/Box
2804SGA-24	Unifusor II with aneroid gauge & stopcock valve	1000mL	24/Box
2805SGA-10	Unifusor II with aneroid gauge & stopcock valve	3000mL	10/Box



2903TGA-24	Unifusor II with piston gauge & thumbwheel valve	500mL	24/Box
2904TGA-24	Unifusor II with piston gauge & thumbwheel valve	1000mL	24/Box
2905TGA-10	Unifusor II with piston gauge & thumbwheel valve	3000mL	10/Box



2903SGA-24	Unifusor II with piston gauge & stopcock valve	500mL	24/Box
2904SGA-24	Unifusor II with piston gauge & stopcock valve	1000mL	24/Box
2905SGA-10	Unifusor II with piston gauge & stopcock valve	3000mL	10/Box

Disposable Replacement Parts




2803X-24	Unifusor II disposable white infusion bag	500mL	24/Box
2804X-24	Unifusor II disposable white infusion bag	1000mL	24/Box
2805X-10	Unifusor II disposable white infusion bag	500mL	10/Box

Reusable Replacement Parts



899F-04	Aneroid gauge & thumbwheel valve assembly	4/Pk
899S-04	Aneroid gauge & stopcock valve assembly	4/Pk



999T-04	Piston gauge & thumbwheel valve assembly	4/Pk
999S-04	Piston gauge & stopcock valve assembly	4/Pk

Unifusor® Pressure Infusors

Unifusor Pressure Infusors are available in 3 Models: Unifusor Disposable, Unifusor II & Unifusor Plus Reusable. All models are available in 500, 100 and 3,000 mL sizes. The Unifusors are 100% tested for leakage during the manufacturing process. Unifusor pressure Infusors include an IV pole loop, a fluid bag hood and a reinforced translucent sleeve which enables easy viewing of fluid type and level. This unit is designed for rapid IV infusion, invasive pressure monitoring and auto-transfusion procedures.



Order Code	Description	Size	Qty
------------	-------------	------	-----

Unifusor Plus Reusable Pressure Infusor



Reusable Infusor - blue bag with white hose

1103FGA-24	Unifusor Plus with aneroid gauge & thumbwheel valve	500mL	24/Box
1104FGA-24	Unifusor Plus with aneroid gauge & thumbwheel valve	1000mL	24/Box
1105FGA-10	Unifusor Plus with aneroid gauge & thumbwheel valve	3000mL	10/Box



1103SGA-24	Unifusor Plus with aneroid gauge & stopcock valve	500mL	24/Box
1104SGA-24	Unifusor Plus with aneroid gauge & stopcock valve	1000mL	24/Box
1105SGA-10	Unifusor Plus with aneroid gauge & stopcock valve	3000mL	10/Box



1103TGA-24	Unifusor Plus with piston gauge & thumbwheel valve	500mL	24/Box
1104TGA-24	Unifusor Plus with piston gauge & thumbwheel valve	1000mL	24/Box
1105TGA-10	Unifusor Plus with piston gauge & thumbwheel valve	3000mL	10/Box



1113SGA-24	Unifusor Plus with piston gauge & stopcock valve	500mL	24/Box
1114SGA-24	Unifusor Plus with piston gauge & stopcock valve	1000mL	24/Box
1115SGA-10	Unifusor Plus with piston gauge & stopcock valve	3000mL	10/Box

Reusable replacement parts



1103XB-05	Unifusor Plus blue reusable bag with white hose	500mL	5/Box
1104XB-05	Unifusor Plus blue reusable bag with white hose	1000mL	5/Box
1105XB-05	Unifusor Plus blue reusable bag with white hose	3000mL	5/Box



899F-04	Aneroid gauge & thumbwheel valve assembly		4/Pk
899S-04	Aneroid gauge & stopcock valve assembly		4/Pk



999T-04	Piston gauge & thumbwheel valve assembly		4/Pk
999S-04	Piston gauge & stopcock valve assembly		4/Pk

Unifusor® Pressure Infusors

Unifusor Pressure Infusors are available in 3 Models: Unifusor Disposable, Unifusor II & Unifusor Plus Reusable. All models are available in 500, 100 and 3,000 mL sizes. The Unifusors are 100% tested for leakage during the manufacturing process. Unifusor pressure Infusors include an IV pole loop, a fluid bag hood and a reinforced translucent sleeve which enables easy viewing of fluid type and level. This unit is designed for rapid IV infusion, invasive pressure monitoring and auto-transfusion procedures.



Order Code	Description	Size	Qty
------------	-------------	------	-----

Unifusor Plus Reusable Pressure Infusor



Reusable Infusor - blue bag with black hose

1103	Unifusor Plus with trigger / gauge assembly	500mL	1 Each
1104	Unifusor Plus with trigger / gauge assembly	1000mL	1 Each
1105	Unifusor Plus with trigger / gauge assembly	3000mL	1 Each

Reusable replacement parts



1103X	Unifusor Plus blue reusable bag with black hose	500mL	1 Each
1104X	Unifusor Plus blue reusable bag with black hose	1000mL	1 Each
1105X	Unifusor Plus blue reusable bag with black hose	3000mL	1 Each



1199	Trigger/ gauge assembly only		1 Each
------	------------------------------	--	--------