




SUPPLIES & ACCESSORIES PULSE OXIMETRY



TruLink SpO2 Supplies & Accessories
Spacelabs TruLink Technology

TruLink branded pulse oximetry supplies and accessories are available for Spacelabs command modules (90496-U/91496-U). These cables and sensors are designed specifically for the Spacelabs module and provide a high level of validated reliable and accurate data capture. Economical and high quality, TruLink SPO2 supplies are available in reusable and disposable sensors for all patient types, including a popular reusable multi-site sensor that depending on patient type can read from finger, ear, foot or toe.











Order Code	Description	Length/Size	Qty
Spacelabs TruLink Technology			
Cables			
 700-0030-00	TruLink SpO ₂ cable	305 cm/10 ft.	1 Each
Reusable sensors			
 015-0660-00	TruLink SpO ₂ finger sensor	Adult >30 kg/66 lbs.	1 Each
 015-0661-00	TruLink SpO ₂ Y multi-site sensor <u>Includes:</u> 1 Sensor 1 ear clip 10 adult finger/neonate feet/hand wraps 10 pediatric finger or infant toe adhesive wraps	Universal Adult >30 kg Pediatric 10-50kg Infant 3-20kg Neonatal <3kg	1 Each

TruLink SpO₂ Supplies & Accessories
Spacelabs TruLink Technology

TruLink branded pulse oximetry supplies and accessories are available for Spacelabs command modules (90496-U/91496-U). These cables and sensors are designed specifically for the Spacelabs module and provide a high level of validated reliable and accurate data capture. Economical and high quality, TruLink SPO₂ supplies are available in reusable and disposable sensors for all patient types, including a popular reusable multi-site sensor that depending on patient type can read from finger, ear, foot or toe.



Order Code	Description	Length/Size	Qty
Spacelabs TruLink Technology			
Disposable Sensors			
	015-0662-00 TruLink SpO ₂ foam sensor	Adult > 30 kg/66 lbs.	24/Box
	015-0663-00 TruLink SpO ₂ vinyl sensor	Adult >30 kg/66 lbs.	24/Box
	015-0664-00 TruLink SpO ₂ foam sensor	Pediatric 10-50 kg/ 22-110 lbs.	24/Box
	015-0665-00 TruLink SpO ₂ vinyl sensor	Pediatric 10-50 kg/ 22-110 lbs.	24/Box
	015-0666-00 TruLink SpO ₂ cloth sensor	Infant 3-20 kg/7-44 lbs.	24/Box
	015-0667-00 TruLink SpO ₂ cloth sensor	Neonatal <3 kg/7 lbs.	24/Box
Dixtal Cables			
	700-0029-00 TruLink SpO ₂ cable compatible with Novametrix® sensors	30 cm/12 in.	1 Each
Dixtal Reusable Sensors			
	015-0321-01 Dixtal SpO ₂ universal sensor ear clip	Universal >30 kg/66 lbs. or <30 kg/66 lbs.	1 Each
Accessories			
	015-0330-00 SpO ₂ foam wrap	Universal	25/Pack

Masimo® SpO2 Supplies & Accessories
LNCS® Masimo SET® Technology

For use with Ultraview SL™ Command Module (91496-M).
 Cables and connectors are color-coded gray.



Order Code	Description	Length/Size	Qty
LNCS Masimo SET Technology			
LNCS Adapter Cable			
700-0906-01	TruLink SpO ₂ cable compatible with Masimo LNCS sensors (91496-M/91496-U)	305 cm/10 ft	1 Each
LNCS Reusable Sensors			
690-0247-00	Masimo LNCS DC-I, SpO ₂ digit sensor	Adult > 30 kg	1 Each
690-0248-00	Masimo LNCS DC-IP, SpO ₂ digit sensor	Pediatric 10-50 kg	1 Each
690-0249-00	Masimo LNCS TC-I, Tip-Clip SpO ₂ ear sensor	Adult > 30 kg	1 Each
690-0250-00	Masimo LNCS TF-I, Transflectance SpO ₂ forehead sensor	Adult > 30 kg	1 Each

Masimo® SpO2 Supplies & Accessories
LNCS® Masimo SET® Technology

For use with Ultraview SL™ Command Module (91496-M).
 Cables and connectors are color-coded gray.



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

LNCS Masimo SET Technology

LNCS Disposable Sensors



690-0251-00	Masimo LNCS Aadx, adhesive SpO2sensor, 18 in.	Adult > 30 kg	20/Box
690-0252-00	Masimo LCNS Pdx, adhesive SpO2 sensor, 18 in.	Pediatric 10 - 50 kg	20/Box
690-0258-00	Masimo LNCS Inf-L, adhesive SpO2 sensor, 3 ft.	Infant 3 kg > 20 kg	20/Box
690-0259-00	Masimo LNCS Neo-L, adhesive SpO2 sensor, 3 ft.	Neonatal < 3 kg or > 40 kg	20/Box

Masimo® SpO₂ Supplies & Accessories
RD SET Technology

For use with Spacelabs Command Module (91496-M) .
 Cables and connectors are color-coded gray.



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

RD SET Masimo Technology



RD SET Adapter Cable

010-1649-00	Masimo RD SET SpO ₂ Cable	12 ft.	1 Each
-------------	--------------------------------------	--------	--------



RD SET Reusable Sensors

690-0267-00	Masimo RD SET DCI reusable SpO ₂ sensor	Adult > 30 kg	1 Each
690-0270-00	Masimo RD SET DCI-P reusable SpO ₂ sensor	Pediatric 10-50 kg	1 Each
690-0274-00	Masimo RD SET TC-I Tip-Clip reusable SpO ₂ sensor	Adult > 30 kg	20/Box

Masimo® SpO₂ Supplies & Accessories
RD SET Technology






For use with Spacelabs Command Module (91496-M) .
 Cables and connectors are color-coded gray.



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

RD SET Masimo Technology

RD SET Adhesive Disposable Sensors

	690-0283-00	Masimo RD SET Adt, adhesive SpO ₂ sensor	Adult	20/Box
	690-0284-00	Masimo RD SET Pdt, adhesive SpO ₂ sensor	Pediatric	20/Box
	690-0285-00	Masimo RD SET Inf, adhesive SpO ₂ sensor	Infant	20/Box
	690-0286-00	Masimo RD SET Neo, adhesive SpO ₂ sensor	Neo/Adult	20/Box
	690-0287-00	Masimo RD SET NeoPt, adhesive SpO ₂ sensor	Neonatal	20/Box

Nellcor® SpO₂ Supplies & Accessories
Nellcor™ OxiMax® Technology

For use with Spacelabs Command Modules (90496,91496-N & 91496-U) and mCare Monitor 91220 and DM3 Monitor.
 Cables and connectors are color-coded purple.



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

Nellcor OxiMax Technology

Cable



700-0792-00	TruLink SpO ₂ cable compatible with Nellcor OxiMax sensors (91496-N, 91220, 91330-N)		1 Each
-------------	---	--	--------



700-0030-00	TruLink SpO ₂ cable compatible with Nellcor OxiSensor II and OxiMax sensors (90496, 91496-U, 91330-S)		1 Each
-------------	--	--	--------

Reusable Sensors

Note: Nellcor sensors are available for purchase through Spacelabs Healthcare only with the purchase of capital equipment. All future purchases should be made directly through Nellcor(Covidien), now part of Medtronic.



690-0003-01	Nellcor DS-100A, finger SpO ₂ sensor	Adult	1 Each
-------------	---	-------	--------



690-0004-01	Nellcor OXI-A/N, SpO ₂ sensor	Adult/Neonatal	1 Each
-------------	--	----------------	--------



690-0039-01	Nellcor OXI-P/I, SpO ₂ sensor	Pediatric/Infant	1 Each
-------------	--	------------------	--------



690-0215-01	Nellcor DS-Y, multi-site SpO ₂ sensor	Universal	1 Each
-------------	--	-----------	--------