

Vital Sign Monitoring

Ultraview DM4 supports measurement of SpO2, pulse rate, NIBP and temperature, and can be used in any patient care environment where basic vital signs monitoring is needed.



91331

Order Code	Description	Size	Length/Size	Qty
------------	-------------	------	-------------	-----

TruLink Nylon Reusable BP Cuffs



Single Tube				
715-1121-10	Single tube, Bayonet	Child	13 - 20 cm	1 Each
715-1131-10	Single tube, Bayonet	Small Adult	18 - 26 cm	1 Each
715-1141-10	Single tube, Bayonet	Adult	26 - 35 cm	1 Each
715-1151-10	Single tube, Bayonet	Adult Long	29 - 38 cm	1 Each
715-1161-10	Single tube, Bayonet	Large Adult	32 - 42 cm	1 Each
715-1171-10	Single tube, Bayonet	Large Adult Long	35 - 44 cm	1 Each
715-1181-10	Single tube, Bayonet	Thigh Cuff	42 - 50 cm	1 Each

TruLink Soft Disposable BP Cuffs



Single Tube				
715-1123-10	Single tube, Bayonet	Child	13 - 20 cm	5/Bag
715-1133-10	Single tube, Bayonet	Small Adult	18 - 26 cm	5/Bag
715-1143-10	Single tube, Bayonet	Adult	26 - 35 cm	5/Bag
715-1153-10	Single tube, Bayonet	Adult Long	29 - 38 cm	5/Bag
715-1163-10	Single tube, Bayonet	Large Adult	32 - 42 cm	5/Bag
715-1173-10	Single tube, Bayonet	Large Adult Long	35 - 44 cm	5/Bag
715-1183-10	Single tube, Bayonet	Thigh Cuff	42 - 50 cm	5/Bag

TruLink Vinyl Disposable BP Cuffs



Single Tube				
715-1122-10	Single tube, Bayonet	Child	13 - 20 cm	5/Bag
715-1132-10	Single tube, Bayonet	Small Adult	18 - 26 cm	5/Bag
715-1142-10	Single tube, Bayonet	Adult	26 - 35 cm	5/Bag
715-1152-10	Single tube, Bayonet	Adult Long	29 - 38 cm	5/Bag
715-1162-10	Single tube, Bayonet	Large Adult	32 - 42 cm	5/Bag
715-1172-10	Single tube, Bayonet	Large Adult Long	35 - 44 cm	5/Bag
715-1182-10	Single tube, Bayonet	Thigh Cuff	42 - 50 cm	5/Bag



Single Tube Hose				
714-0018-00	Single tube, quick connect	Adult	275 cm/9 ft	1 Each
714-0060-00	Single tube, quick connect	Adult	366 cm/12 ft	1 Each

Vital Sign Monitoring

Ultraview DM4 supports measurement of SpO2, pulse rate, NIBP and temperature, and can be used in any patient care environment where basic vital signs monitoring is needed.



91331

Order Code	Description	Size	Length/Size	Qty
------------	-------------	------	-------------	-----

SoftCheck Fabric Neonatal Disposable BP Cuffs



Plastic Bayonet

Single Tube					
NN1ST-HP-10	Single tube, HP	Size 1	3-6 cm	10/Box	
NN2ST-HP-10	Single tube, HP	Size 2	4-8 cm	10/Box	
NN3ST-HP-10	Single tube, HP	Size 3	6-11 cm	10/Box	
NN4ST-HP-10	Single tube, HP	Size 4	7-13 cm	10/Box	
NN5ST-HP-10	Single tube, HP	Size 5	8-15 cm	10/Box	

SoftCheck Vinyl Neonatal Disposable BP Cuffs



Plastic Bayonet

Single Tube					
VNN1ST-HP-10	Single tube, HP	Size 1	3-6 cm	10/Box	
VNN2ST-HP-10	Single tube, HP	Size 2	4-8 cm	10/Box	
VNN3ST-HP-10	Single tube, HP	Size 3	6-11 cm	10/Box	
VNN4ST-HP-10	Single tube, HP	Size 4	7-13 cm	10/Box	
VNN5ST-HP-10	Single tube, HP	Size 5	8-15 cm	10/Box	

Vital Sign Monitoring

Ultraview DM4 supports measurement of SpO₂, pulse rate, NIBP and temperature, and can be used in any patient care environment where basic vital signs monitoring is needed.



91331

Order Code	Description	Size	Length/Size	Qty
------------	-------------	------	-------------	-----

LNCS Masimo SET SpO₂ Technology - 91331 with SPO2 Option M

LNCS Adapter Cable

010-2227-00	DM4 Masimo LNCS SpO ₂ cable compatible with Masimo LNCS sensors		305 cm/10 ft	1 Each
-------------	--	--	--------------	--------

LNCS reusable sensors



690-0247-00	Masimo LNCS DC-I Digit Sensor		Adult > 30 kg	1 Each
690-0248-00	Masimo LNCS DC-IP Pediatric Digit Sensor		Pediatric 10-50 kg	1 Each
690-0249-00	Masimo LNCS TC-I Tip-Clip Ear Sensor		Adult > 30 kg	1 Each
690-0250-00	Masimo LNCS TF-I Transflectance Forehead Sensor		Adult > 30 kg	1 Each
690-0291-00	Masimo LNCS Y-I multisite sensor		> 1 kg	1 Each

Vital Sign Monitoring

Ultraview DM4 supports measurement of SpO2, pulse rate, NIBP and temperature, and can be used in any patient care environment where basic vital signs monitoring is needed.



91331

Order Code	Description	Size	Length/Size	Qty
------------	-------------	------	-------------	-----

LNCS Masimo SET SpO2 Technology - 91331 with SPO2 Option M

LNCS Disposable Sensors



690-0251-00	Masimo LNCS Adtx, adhesive SpO2sensor, 18 in.	Adult > 30 kg	20/Box
690-0252-00	Masimo LCNS Pdtx, adhesive SpO2 sensor, 18 in.	Pediatric 10 - 50 kg	20/Box
690-0258-00	Masimo LNCS Inf-L, adhesive SpO2 sensor, 3 ft.	Infant 3 kg > 20 kg	20/Box
690-0259-00	Masimo LNCS Neo-L, adhesive SpO2 sensor, 3 ft.	Neonatal < 3 kg or > 40 kg	20/Box

Vital Sign Monitoring

Ultraview DM4 supports measurement of SpO2, pulse rate, NIBP and temperature, and can be used in any patient care environment where basic vital signs monitoring is needed.



91331

Order Code	Description	Size	Length/Size	Qty
------------	-------------	------	-------------	-----

Masimo RD SET SpO2 Technology - 91331 with SPO2 Option M

RD SET Adapter Cable



010-2193-00	DM4 Masimo RD SET SpO2 cable compatible with Masimo RD SET sensors		305 cm/10 ft	1 Each
-------------	--	--	--------------	--------

RD SET Reusable Sensors








690-0267-00	Masimo RD SET DCI reusable SpO2 sensor		Adult > 30 kg	1 Each
690-0270-00	Masimo RD SET DCI-P reusable SpO2 sensor		Pediatric 10-50 kg	1 Each
690-0274-00	Masimo RD SET TC-I Tip-Clip reusable SpO2 sensor		Adult > 30 kg	1 Each
690-0281-00	Masimo RD SET Y-I multi-site reusable SpO2 sensor		> 1 kg	1 Each

Vital Sign Monitoring

Ultraview DM4 supports measurement of SpO2, pulse rate, NIBP and temperature, and can be used in any patient care environment where basic vital signs monitoring is needed.



91331

Order Code	Description	Size	Length/Size	Qty
Masimo RD SET SpO2 Technology - 91331 with SPO2 Option M				
RD SET Adhesive Disposable Sensors				
 RD SET Adt	690-0283-00	Masimo RD SET Adt, SpO2 adhesive Sensor	Adult > 30 kg finger or toe	20/Box
 RD SET Pdt	690-0284-00	Masimo RD SET Pdt, SpO2 adhesive Sensor	Pediatric 10-50 kg finger or toe	20/Box
 RD SET Inf	690-0285-00	Masimo RD SET Inf, SpO2 adhesive Sensor	Infant 3—20 kg thumb or great toe	20/Box
 RD SET Neo	690-0286-00	Masimo RD SET Neo, SpO2 adhesive Sensor	Neonatal <3 kg, hand or foot Adult > 40 kg finger or toe	20/Box
 RD SET NeoPt	690-0287-00	Masimo RD SET NeoPt, SpO2 adhesive Sensor	Neonatal <1 kg hand or foot	20/Box

Vital Sign Monitoring

Ultraview DM4 supports measurement of SpO2, pulse rate, NIBP and temperature, and can be used in any patient care environment where basic vital signs monitoring is needed.



91331

Order Code	Description	Size	Length/Size	Qty
Nellcor OxiMax SpO2 Technology - 91331 with SPO2 Option N				
Cable				
010-2192-00	DM4 Nellcor OxiMax SpO2 cable compatible with Nellcor OxiMax sensors			1 Each
Reusable Sensors				
Note: Nellcor sensors are available for purchase through Spacelabs Healthcare only with the purchase of capital equipment. All future purchases should be made directly through Nellcor(Covidien).				
690-0003-01	Nellcor OxiMax DS-100A finger sensor		Adult	1 Each
690-0004-01	Nellcor OXI-A/N sensor		Adult/Neonatal	1 Each
690-0039-01	Nellcor OXI-P/I sensor		Pediatric/Infant	1 Each
690-0215-01	Nellcor DS-Y multi-site sensor		Universal	1 Each

Vital Sign Monitoring

Ultraview DM4 supports measurement of SpO2, pulse rate, NIBP and temperature, and can be used in any patient care environment where basic vital signs monitoring is needed.



91331

Order Code	Description	Size	Length/Size	Qty
------------	-------------	------	-------------	-----

Temperature

Filac



010-2175-00	Filac temperature probe			1 Each
-------------	-------------------------	--	--	--------



010-2174-00	Filac temperature isolation chamber			1 Each
-------------	-------------------------------------	--	--	--------



010-2228-00	Filac temperature oral probe covers, 20/bx			25/bxs
-------------	--	--	--	--------

Exergen



010-2157-00	Temporal Artery Thermometer			1 Each
-------------	-----------------------------	--	--	--------