


Temperature Supplies

A variety of supplies are available through Spacelabs for 90496/91496 modules configured for temperature. A selection of skin, tympanic, esophageal and rectal temperature probes are available in different diameters and lengths.



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

Continuous Temperature

Cables

For use with Command Modules 90496/91496, DM3 Monitor 91330 and mCare Monitor 91220



700-0031-00	TruLink Dual Temperature Cable	30 cm/12 in.	1 Each
-------------	--------------------------------	--------------	--------



700-4000-00	Universal temperature probe extension cable, YSI 400 series compatible *Required for use of single use temperature probes/sensors	305 cm/10 ft	1 Each
-------------	--	--------------	--------

Reusable Probes/Sensors

Skin Surface



690-0028-01	YSI 700 series compatible with 6.5 mm DIA phone jack type connections	305 cm/10 ft	1 Each
-------------	---	--------------	--------

Esophageal / Rectal




690-0272-00	12 French, YSI 400 series compatible	305 cm/10 ft	1 Each
-------------	--------------------------------------	--------------	--------









118-0256-01	12 French, YSI 700 series compatible	305 cm/10 ft	1 Each
-------------	--------------------------------------	--------------	--------

Temperature Supplies

A variety of supplies are available through Spacelabs for 90496/91496 modules configured for temperature. A selection of skin, tympanic, esophageal and rectal temperature probes are available in different diameters and lengths.



Order Code	Description	Length/Size	Qty
Single use Probes/Sensors - Requires YSI 400 series cable (700-4000-00)			
Skin			
 701-4000-00	Skin Temperature Sensor		20/Box
Esophageal/ Stethoscope			
 701-4009-00	9 French, YSI 400 series compatible		20/Box
701-4012-00	12 French, YSI 400 series compatible		20/Box
701-4018-00	18 French, YSI 400 series compatible		20/Box
Esophageal/ Rectal			
 701-4109-00	9 French, YSI 400 series compatible		20/Box
701-4112-00	12 French, YSI 400 series compatible		20/Box
Foley catheter with temperature probe			
 701-4214-00	14 French, YSI 400 series compatible		20/Box
701-4216-00	16 French, YSI 400 series compatible		20/Box
701-4218-00	18 French, YSI 400 series compatible		20/Box
Tympanic temperature sensor			
 701-4300-00	Adult, YSI 400 series compatible		20/Box
Exergen			
<i>For use with Qube Mini (91389), Qube (91390), Xprezzon(91393) & DM4 (91331) Monitors</i>			
 010-2157-00	Temporal Artery Thermometer		1 Each